

SKU #	Product Name	Negative* Calibrator	Low Calibrator	Level 1* Control	Cutoff* Calibrator	Level 2* Control	Intermed. Cal #1	Intermed. Cal #2	High Cal
60-AC0580-100 60-AC0580-1L	6-Acetylmorphine/ Heroin Metabolite	0 (60-AC0598-5)	5 (60-AC0584-5)	7.5 (60-AC0594-5)	10 (60-AC0586-5)	12.5 (60-AC0596-5)	20 (60-AC0588-5)		40 (60-AC0590-5)
60-AM0200-100 60-AM0200-1L	Amphetamines 500 (d-Methamphetamine)	0 (60-UN0002-5)	250 (60-AM0204-5)	375 (60-AM0214-5)	500 (60-AM0206-5)	625 (60-AM0216-5)	1000 (60-AM0208-5)		2000 (60-AM0210-5)
60-AM0080-100 60-AM0080-1L	Amphetamines 1000 (d-Methamphetamine)	0 (60-UN0002-5)	500 (60-AM0084-5)	750 (60-AM0094-5)	1000 (60-AM0086-5)	1250 (60-AM0096-5)	1500 (60-AM0088-5)		2000 (60-AM0090-5)
60-BA0280-100 60-BA0280-1L	Barbiturate 200 (Secobarbital)	0 (60-UN0002-5)	100 (60-BA0284-5)	100 (60-BA0284-5)	200 (60-BA0286-5)	300 (60-BA0288-5)	300 (60-BA0288-5)		1000 (60-BA0290-5)
60-BA0280-100 60-BA0280-1L	Barbiturate 300 (Secobarbital)	0 (60-UN0002-5)	100 (60-BA0284-5)	200 (60-BA0286-5)	300 (60-BA0288-5)	400 (60-BA0292-5)	200 (60-BA0286-5)		1000 (60-BA0290-5)
60-BZ0260-100 60-BZ0260-1L	Benzodiazepine 200 (Oxazepam)	0 (60-UN0002-5)	100 (60-BZ0264-5)	100 (60-BZ0264-5)	200 (60-BZ0266-5)	300 (60-BZ0268-5)	300 (60-BZ0268-5)		1000 (60-BZ0270-5)
60-BZ0260-100 60-BZ0260-1L	Benzodiazepine 300 (Oxazepam)	0 (60-UN0002-5)	100 (60-BZ0264-5)	200 (60-BZ0266-5)	300 (60-BZ0268-5)	400 (60-BZ0272-5)	200 (60-BZ0266-5)		1000 (60-BZ0270-5)
60-BP0540-100 60-BP0540-1L	Buprenorphine 5 (Norbuprenorphine)	0 (60-UN0002-5)		3 (60-NB0544-5)	5 (60-NB0546-5)	7 (60-NB0548-5)	20 (60-NB0554-5)	40 (60-NB0556-5)	75 (60-NB0558-5)
60-BP0540-100 60-BP0540-1L	Buprenorphine 10 (Norbuprenorphine)	0 (60-UN0002-5)	5 (60-NB0546-5)	7 (60-NB0548-5)	10 (60-NB0550-5)	13 (60-NB0552-5)	20 (60-NB0554-5)	40 (60-NB0556-5)	75 (60-NB0558-5)
60-CC0680-100 60-CC0680-1L	Cocaine Metabolite 150 (Benzoylcgonine)	0 (60-UN0002-5)	75 (60-CC0684-5)	112.5 (60-CC0694-5)	150 (60-CC0686-5)	187.5 (60-CC0696-5)	300 (60-CC0688-5)		1000 (60-CC0690-5)
60-CC0060-100 60-CC0060-1L	Cocaine Metabolite 300 (Benzoylcgonine)	0 (60-UN0002-5)	150 (60-CC0064-5)	225 (60-CC0074-5)	300 (60-CC0066-5)	375 (60-CC0076-5)	1000 (60-CC0068-5)		3000 (60-CC0070-5)
60-EC0320-100 60-EC0320-1L	Ecstasy/MDMA	0 (60-UN0002-5)	100 (60-EC0324-5)	375 (60-EC0334-5)	500 (60-EC0326-5)	625 (60-EC0336-5)	750 (60-EC0328-5)		1000 (60-EC0330-5)
60-EA0440-100 60-EA0440-1L	Ethyl Alcohol **Quantitative testing only	0 (60-UN0002-5)		50 (60-EA0448-5)	100 (60-EA0446-5)	300 (60-EA0450-5)			
60-ME0220-100 60-ME0220-1L	Methadone 300	0 (60-UN0002-5)	150 (60-ME0224-5)	225 (60-ME0234-5)	300 (60-ME0226-5)	375 (60-ME0236-5)	600 (60-ME0228-5)		1000 (60-ME0230-5)
60-MM0380-100 60-MM0380-1L	Methadone Metabolite/EDDP 100	0 (60-UN0002-5)	50 (60-MM0404-5)	75 (60-MM0414-5)	100 (60-MM0406-5)	125 (60-MM0416-5)	250 (60-MM0408-5)		500 (60-MM0410-5)
60-MM0380-100 60-MM0380-1L	Methadone Metabolite/EDDP 300	0 (60-UN0002-5)	150 (60-MM0384-5)	225 (60-MM0394-5)	300 (60-MM0386-5)	375 (60-MM0396-5)	600 (60-MM0388-5)		1000 (60-MM0390-5)
60-MP0700-100 60-MP0700-1L	Methamphetamine (d-Methamphetamine)	0 (60-UN0002-5)	250 (60-MP0704-5)	375 (60-MP0714-5)	500 (60-MP0706-5)	625 (60-MP0716-5)	1000 (60-MP0708-5)		2000 (60-MP0710-5)
60-OP0040-100 60-OP0040-1L	Opiates 300 (Morphine)	0 (60-UN0002-5)	150 (60-OP0044-5)	225 (60-OP0054-5)	300 (60-OP0046-5)	375 (60-OP0056-5)	600 (60-OP0048-5)		1000 (60-OP0050-5)
60-OP0660-100 60-OP0660-1L	Opiates 2000 (Morphine)	0 (60-UN0002-5)	1000 (60-OP0664-5)	1500 (60-OP0674-5)	2000 (60-OP0666-5)	2500 (60-OP0676-5)	4000 (60-OP0668-5)		6000 (60-OP0670-5)
60-OX0480-100 60-OX0480-1L	Oxycodone 100	0 (60-UN0002-5)	50 (60-OX0504-5)	75 (60-OX0484-5)	100 (60-OX0486-5)	125 (60-OX0488-5)	300 (60-OX0492-5)	500 (60-OX0496-5)	800 (60-OX0498-5)
60-OX0480-100 60-OX0480-1L	Oxycodone 300	0 (60-UN0002-5)	50 (60-OX0504-5)	225 (60-OX0490-5)	300 (60-OX0492-5)	375 (60-OX0494-5)	500 (60-OX0496-5)		800 (60-OX0498-5)
60-PH0020-100 60-PH0020-1L	Phencyclidine	0 (60-UN0002-5)	12.5 (60-PC0024-5)	18 (60-PC0034-5)	25 (60-PC0026-5)	32 (60-PC0036-5)	50 (60-PC0028-5)		100 (60-PC0030-5)
60-PP0240-100 60-PP0240-1L	Propoxyphene	0 (60-UN0002-5)	150 (60-PP0244-5)	225 (60-PP0254-5)	300 (60-PP0246-5)	375 (60-PP0256-5)	600 (60-PP0248-5)		1000 (60-PP0250-5)
60-CA0140-100 60-CA0140-1L	THC/Cannabinoid 25	0 (60-TH0004-5)	12.5 (60-TH0144-5)	18.75 (60-TH0012-5)	25 (60-TH0146-5)	31.25 (60-TH0148-5)	37.5 (60-TH0014-5)		50 (60-TH0150-5)
60-CA0140-100 60-CA0140-1L	THC/Cannabinoid 50	0 (60-TH0004-5)	25 (60-TH0146-5)	37.5 (60-TH0014-5)	50 (60-TH0150-5)	62.5 (60-TH0016-5)	75 (60-TH0152-5)		100 (60-TH0154-5)
60-CA0140-100 60-CA0140-1L	THC/Cannabinoid 100	0 (60-TH0004-5)	50 (60-TH0150-5)	75 (60-TH0152-5)	100 (60-TH0154-5)	125 (60-TH0018-5)	150 (60-TH0156-5)		200 (60-TH0158-5)
60-HY0760-100 60-HY0760-1L	Hydrocodone 100	0 (60-UN0002-5)	50 (60-HY0764-5)	75 (60-HY0774-5)	100 (60-HY0764-5)	125 (60-HY0776-5)	150 (60-HY0768-5)		300 (60-HY0770-5)
60-HY0780-100 60-HY0780-1L	Hydrocodone 300	0 (60-UN0002-5)	150 (60-HY0784-5)	225 (60-HY0794-5)	300 (60-HY0786-5)	375 (60-HY0796-5)	500 (60-HY0788-5)		800 (60-HY0790-5)
60-FN0620-100 60-FN0620-1L	Fentanyl (Norfentanyl)	0 (60-UN0002-5)		3.75 (60-NF0634-5)	5 (60-NF0626-5)	6.25 (60-NF0636-5)			
60-CM0740-100 60-CM0740-1L	Carisoprodol Metabolite (Meprobamate)	0 (60-UN0002-5)	50 (60-MB0744-5)	75 (60-MB0754-5)	100 (60-MB0746-5)	125 (60-MB0756-5)	200 (60-MB0748-5)		400 (60-MB0750-5)

* Use this combination of Calibrators and Controls for Qualitative Tests

SKU #	Multi-Analyte "C" with 6 Analytes	Universal Neg. Cal. 60-UN0002-5	Low Calibrator 60-MA1704-15	Level 1 Control 60-MA1714-15	Cutoff Calibrator 60-MA1706-15	Level 2 Control 60-MA1716-15	Intermed. Calibrator 60-MA1708-15	High Calibrator 60-MA1710-15
60-CC0680	Cocaine Metabolite 150 (Benzoyllecgonine)	0 ng/mL	75 ng/mL	112.5 ng/mL	150 ng/mL	187.5 ng/mL	300 ng/mL	1000 ng/mL
60-AM0200 60-MP0700	Amphetamines 500 & Methamphetamine 500 (d-Methamphetamine)	0 ng/mL	250 ng/mL	375 ng/mL	500 ng/mL	625 ng/mL	750 ng/mL	1000 ng/mL
60-ME0220	Methadone 300	0 ng/mL	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	600 ng/mL	1000 ng/mL
60-OP0660 60-OP0660	Opiates 2000 (Morphine)	0 ng/mL	1000 ng/mL	1500 ng/mL	2000 ng/mL	2500 ng/mL	4000 ng/mL	6000 ng/mL
60-BZ0260	Benzodiazepine 200 (Oxazepam)	0 ng/mL	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL
60-BA0280	Barbiturate 200 (Secobarbital)	0 ng/mL	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL

SKU #	Multi-Analyte "D" with 5 Analytes	Negative Calibrator 60-MA1742-15	Low Calibrator 60-MA1744-15	Level 1 Control 60-MA1754-15	Cutoff Calibrator 60-MA1746-15	Level 2 Control 60-MA1756-15	Intermed. Calibrator 60-MA1748-15	High Calibrator 60-MA1750-15
60-CC0060	Cocaine Metabolite 300 (Benzoyllecgonine)	0 ng/mL	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	1000 ng/mL	3000 ng/mL
60-AM0200 60-AM0200	Amphetamines 1000	0 ng/mL	500 ng/mL	750 ng/mL	1000 ng/mL	1250 ng/mL	1500 ng/mL	2000 ng/mL
60-OP0040	Opiates 300 (Morphine)	0 ng/mL	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	600 ng/mL	1000 ng/mL
60-BZ0260	Benzodiazepine 200 (Oxazepam)	0 ng/mL	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL
60-OX0480	Oxycodone 100**	0 ng/mL	50 ng/mL	75 ng/mL	100 ng/mL	125 ng/mL	300 ng/mL	500 ng/mL

** Requires 800 ng/mL high calibrator purchased separately

SKU #	Multi-Analyte "E" with 4 Analytes	Negative Calibrator 60-MA1782-15	Low Calibrator 60-MA1784-15	Level 1 Control 60-MA1794-15	Cutoff Calibrator 60-MA1786-15	Level 2 Control 60-MA1796-15	Intermed. Calibrator 60-MA1788-15	High Calibrator 60-MA1790-15
60-BA0280	Barbiturate 200 (Secobarbital)	0 ng/mL	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL
60-BP0540	Buprenorphine 10*** (Norbuprenorphine)	0 ng/mL	5 ng/mL	7 ng/mL	10 ng/mL	13 ng/mL	20 ng/mL	40 ng/mL
60-ME0220	Methadone 300	0 ng/mL	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	600 ng/mL	1000 ng/mL
60-PH0020	Phencyclidine 25	0 ng/mL	12.5 ng/mL	18 ng/mL	25 ng/mL	32 ng/mL	50 ng/mL	100 ng/mL

*** Requires 75 ng/mL high calibrator purchased separately